




2-X Specifications

P/N: IV2SWBK1

Version : 1.0

Performance	
Internal volume	59 μ L
Dead volume	none
Switching time from Port to Port	3 ms
Mechanical	
Weight	36 g
Dimensions	56.4 x 46 x 11 mm
Wetted materials	PEEK, EPDM or FFKM (on demand)
Channel diameter	1.2 mm
Fittings	$\frac{1}{4}$ -28 (1/16" OD) flangeless with flat ferrule
Max pressure	-1 to 2 bars
Operating temperature	0-50 °C (32-122 °F)
Gas compatibility	Resistant to neutral and aggressive gases
Liquid compatibility	Resistant to neutral and aggressive liquids
Electrical	
Power supply	24 VDC
Port communication	RJ45
Software control	OxyGen, SDK